

NF-H Polyclonal Antibody FITC Conjugated

Catalog No: #C90286FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NF-H Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rat
Immunogen Type	1-100/1026
Conjugates	FITC
Target Name	NF-H
Other Names	NF-H
Accession No.	Uniprot:P12036Gene ID:4744
Uniprot	P12036
GeneID	4744
Calculated MW	118
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only