

NF-M Polyclonal Antibody APC Conjugated

Catalog No: #C90288APC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NF-M Polyclonal Antibody APC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rat
Immunogen Type	101-200/916
Conjugates	APC
Target Name	NF-M
Other Names	NF-M
Accession No.	Uniprot:P12839Gene ID:4741
Uniprot	P12839
GeneID	4741
Calculated MW	102
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=2ug/Test IF=1:50-200

Note: This product is for in vitro research use only