

Beta Tubulin Polyclonal Antibody Gold Conjugated

Catalog No: #C90289Gold

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Beta Tubulin Polyclonal Antibody Gold Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IEM IGS GICA
Species Reactivity	Hu Ms Rat
Immunogen Type	31-130/1444
Conjugates	Gold
Target Name	Beta Tubulin
Other Names	Beta Tubulin
Accession No.	Uniprot:Q13509Gene ID:203068
Uniprot	Q13509
GeneID	203068
Calculated MW	50
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IEM=1:20-200 IGS=1:20-200 GICA=1:20-200

Note: This product is for in vitro research use only