

RANKL/CD254 Polyclonal Antibody PE Conjugated

Catalog No: #C90302PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	RANKL/CD254 Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rat
Immunogen Type	210-317/317
Conjugates	PE
Target Name	RANKL/CD254
Other Names	RANKL/CD254
Accession No.	Uniprot:O14788Gene ID:8600
Uniprot	O14788
GeneID	8600
Calculated MW	35
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only