

GAPDH Polyclonal Antibody RBITC Conjugated

Catalog No: #C90304RBITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GAPDH Polyclonal Antibody RBITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	231-335/335
Conjugates	RBITC
Target Name	GAPDH
Other Names	GAPDH
Accession No.	Uniprot:P04406Gene ID:2597
Uniprot	P04406
GeneID	2597
Calculated MW	37
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only