

## Vimentin Polyclonal Antibody Biotin Conjugated

Catalog No: #C90305Bio

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Vimentin Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC-P IHC-F ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	371-466/466
Conjugates	Bio
Target Name	Vimentin
Other Names	Vimentin
Accession No.	Uniprot:P08670Gene ID:7431
Uniprot	P08670
GeneID	7431
Calculated MW	53
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only