

CD4 Polyclonal Antibody PE-Cy7 Conjugated

Catalog No: #C90307PE-Cy7

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD4 Polyclonal Antibody PE-Cy7 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rat
Immunogen Type	231-330/457
Conjugates	PE-Cy7
Target Name	CD4
Other Names	CD4
Accession No.	Uniprot:P06332Gene ID:12504
Uniprot	P06332
GeneID	12504
Calculated MW	48
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=2ug/Test IF=1:50-200

Note: This product is for in vitro research use only