

## CD90/Thy-1 Polyclonal Antibody FITC Conjugated

Catalog No: #C90311FITC

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CD90/Thy-1 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rat
Immunogen Type	31-120/161
Conjugates	FITC
Target Name	CD90/Thy-1
Other Names	CD90/Thy-1
Accession No.	Uniprot:P01830Gene ID:24832
Uniprot	P01830
GeneID	24832
Calculated MW	12
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only