

## IL-8/CXCL8 Polyclonal Antibody HRP Conjugated

Catalog No: #C90312HRP

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	IL-8/CXCL8 Polyclonal Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC-P IHC-F
Species Reactivity	Hu
Immunogen Type	31-99/99
Conjugates	HRP
Target Name	IL-8/CXCL8
Other Names	IL-8/CXCL8
Accession No.	Uniprot:P10145Gene ID:3576
Uniprot	P10145
GeneID	3576
Calculated MW	44417
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IHC-P=1:50-200 IHC-F=1:50-200

Note: This product is for in vitro research use only