

IL-8/CXCL8 Polyclonal Antibody FITC Conjugated

Catalog No: #C90312FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	IL-8/CXCL8 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu
Immunogen Type	31-99/99
Conjugates	FITC
Target Name	IL-8/CXCL8
Other Names	IL-8/CXCL8
Accession No.	Uniprot:P10145Gene ID:3576
Uniprot	P10145
GeneID	3576
Calculated MW	44417
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only