

HIV gp160 Polyclonal Antibody FITC Conjugated

Catalog No: #C90315FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	HIV gp160 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	HIV1
Immunogen Type	661-710/876
Conjugates	FITC
Target Name	HIV gp160
Other Names	HIV gp160
Accession No.	Uniprot:Q79670Gene ID:
Uniprot	Q79670
Calculated MW	95
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only