

## CGRP-I/II Polyclonal Antibody Biotin Conjugated

Catalog No: #C90317Bio

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CGRP-I/II Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC-P IHC-F IF
Species Reactivity	Hu Ms Rat
Immunogen Type	65-130/130
Conjugates	Bio
Target Name	CGRP-I/II
Other Names	CGRP-I/II
Accession No.	Uniprot:P06881/P10092Gene ID:796/797
Uniprot	P06881/P10092
GeneID	796/797
Calculated MW	3.8
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IHC-P=1:50-200 IHC-F=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only