

CD56 Polyclonal Antibody PE Conjugated

Catalog No: #C90322PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | CD56 Polyclonal Antibody PE Conjugated |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | affinity purified by Protein A |
| Applications | Flow-Cyt ICC IF |
| Species Reactivity | Hu Ms Rat |
| Immunogen Type | 621-720/858 |
| Conjugates | PE |
| Target Name | CD56 |
| Other Names | CD56 |
| Accession No. | Uniprot:P13591Gene ID:4684 |
| Uniprot | P13591 |
| GeneID | 4684 |
| Calculated MW | 92 |
| Formulation | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol. |
| Storage | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only