

IL-1 Beta Polyclonal Antibody FITC Conjugated

Catalog No: #C90325FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	IL-1 Beta Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rat
Immunogen Type	101-200/269
Conjugates	FITC
Target Name	IL-1 Beta
Other Names	IL-1 Beta
Accession No.	Uniprot:P01584Gene ID:3553
Uniprot	P01584
GeneID	3553
Calculated MW	17/32
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only