

Lysozyme Polyclonal Antibody PE Conjugated

Catalog No: #C90328PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Lysozyme Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	purified protein
Conjugates	PE
Target Name	Lysozyme
Other Names	Lysozyme
Accession No.	Uniprot:P61626Gene ID:4069
Uniprot	P61626
GeneID	4069
Calculated MW	17
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only