

## laminin Polyclonal Antibody HRP Conjugated

Catalog No: #C90332HRP

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	laminin Polyclonal Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC-P IHC-F
Species Reactivity	Hu Ms Rat
Immunogen Type	901-1000/3075
Conjugates	HRP
Target Name	laminin
Other Names	laminin
Accession No.	Uniprot:P25391 Gene ID:284217
Uniprot	P25391
GeneID	284217
Calculated MW	198
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IHC-P=1:50-200 IHC-F=1:50-200

Note: This product is for in vitro research use only