

Thrombin receptor/ PAR1 Polyclonal Antibody Biotin Conjugated



Catalog No: #C90336Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Thrombin receptor/ PAR1 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC-P IHC-F IF
Species Reactivity	Hu Ms Rat
Immunogen Type	251-350/425
Conjugates	Bio
Target Name	Thrombin receptor/ PAR1
Other Names	Thrombin receptor/ PAR1
Accession No.	Uniprot:P25116Gene ID:2149
Uniprot	P25116
GeneID	2149
Calculated MW	42
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only