

Nanog Polyclonal Antibody APC Conjugated

Catalog No: #C90337APC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Nanog Polyclonal Antibody APC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Ms Rat
Immunogen Type	71-170/305
Conjugates	APC
Target Name	Nanog
Other Names	Nanog
Accession No.	Uniprot:Q80Z64Gene ID:71950
Uniprot	Q80Z64
GeneID	71950
Calculated MW	34
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only