

Myc tag Polyclonal Antibody HRP Conjugated

Catalog No: #C90340HRP

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Myc tag Polyclonal Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC-P IHC-F
Species Reactivity	Hu
Conjugates	HRP
Target Name	Myc tag
Other Names	Myc tag
Accession No.	Uniprot:P01106Gene ID:4609
Uniprot	P01106
GeneID	4609
Calculated MW	49
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:10000-20000 IHC-P=1:50-200 IHC-F=1:50-200

Note: This product is for in vitro research use only