

phospho-AKT (Ser473) Polyclonal Antibody PE Conjugated

Catalog No: #C90343PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	phospho-AKT (Ser473) Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	QF(p-S)YS
Conjugates	PE
Target Name	phospho-AKT (Ser473)
Other Names	phospho-AKT (Ser473)
Accession No.	Uniprot:P47196Gene ID:207
Uniprot	P47196
GeneID	207
Calculated MW	56
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only