

# beta-Amyloid 1-40 (CT) Polyclonal Antibody Bio Conjugated



Catalog No: #C90344Bio

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	beta-Amyloid 1-40 (CT) Polyclonal Antibody Bio Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC-P IHC-F ICC IF WB
Species Reactivity	Hu Ms Rat
Conjugates	Bio
Target Name	beta-Amyloid 1-40 (CT)
Other Names	beta-Amyloid 1-40 (CT)
Accession No.	Uniprot:P05067Gene ID:351
Uniprot	P05067
GeneID	351
Calculated MW	4.3
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200 IF=1:50-200

WB=1:500-2000

Note: This product is for in vitro research use only