

CD144/VE Cadherin Polyclonal Antibody Biotin Conjugated

Catalog No: #C90345Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD144/VE Cadherin Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC-P IHC-F ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	601-700/784
Conjugates	Bio
Target Name	CD144/VE Cadherin
Other Names	CD144/VE Cadherin
Accession No.	Uniprot:P55284Gene ID:12562
Uniprot	P55284
GeneID	12562
Calculated MW	86
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only