

Measles virus fusion protein Polyclonal Antibody Biotin Conjugated

Catalog No: #C90348Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Measles virus fusion protein Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC-F IF
Species Reactivity	Measles virus
Immunogen Type	451-550/550
Conjugates	Bio
Target Name	Measles virus fusion protein
Other Names	Measles virus fusion protein
Accession No.	Uniprot:P69353Gene ID:1489800
Uniprot	P69353
GeneID	1489800
Calculated MW	57
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC-F=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only