

MV hemagglutinin glycoprotein Polyclonal Antibody Biotin Conjugated

Catalog No: #C90349Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MV hemagglutinin glycoprotein Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC-P IHC-F IF
Species Reactivity	Measles virus
Immunogen Type	101-200/617
Conjugates	Bio
Target Name	MV hemagglutinin glycoprotein
Other Names	MV hemagglutinin glycoprotein
Accession No.	Uniprot:P08362Gene ID:
Uniprot	P08362
Calculated MW	68
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC-P=1:50-200 IHC-F=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only