

GFP Polyclonal Antibody Biotin Conjugated

Catalog No: #C90350Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GFP Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC-P IHC-F IF
Species Reactivity	Belt Jellyfish beltjellyfish
Immunogen Type	25-120/238
Conjugates	Bio
Target Name	GFP
Other Names	GFP
Accession No.	Uniprot:P42212Gene ID:107331116
Uniprot	P42212
GeneID	107331116
Calculated MW	27
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only