

## GFP Polyclonal Antibody Gold Conjugated

Catalog No: #C90350Gold

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	GFP Polyclonal Antibody Gold Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IEM IGS GICA
Species Reactivity	Belt Jellyfish beltjellyfish
Immunogen Type	25-120/238
Conjugates	Gold
Target Name	GFP
Other Names	GFP
Accession No.	Uniprot:P42212Gene ID:107331116
Uniprot	P42212
GeneID	107331116
Calculated MW	27
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IEM=1:20-200 IGS=1:20-200 GICA=1:20-200

Note: This product is for in vitro research use only