

## FSH receptor Polyclonal Antibody PE-Cy7 Conjugated

Catalog No: #C90351PE-Cy7

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	FSH receptor Polyclonal Antibody PE-Cy7 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rat
Immunogen Type	201-300/695
Conjugates	PE-Cy7
Target Name	FSH receptor
Other Names	FSH receptor
Accession No.	Uniprot:P23945Gene ID:2492
Uniprot	P23945
GeneID	2492
Calculated MW	78
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only