

VCAM1 Polyclonal Antibody Biotin Conjugated

Catalog No: #C90358Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	VCAM1 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB
Species Reactivity	Hu Ms Rat
Immunogen Type	651-739/739
Conjugates	Bio
Target Name	VCAM1
Other Names	VCAM1
Accession No.	Uniprot:P19320Gene ID:7412
Uniprot	P19320
GeneID	7412
Calculated MW	81
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200

Note: This product is for in vitro research use only