

Mafa Polyclonal Antibody PE-Cy5.5 Conjugated

Catalog No: #C90361PE-Cy5.5

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mafa Polyclonal Antibody PE-Cy5.5 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	265-359/359
Conjugates	PE-Cy5.5
Target Name	Mafa
Other Names	Mafa
Accession No.	Uniprot:Q8CF90Gene ID:378435
Uniprot	Q8CF90
GeneID	378435
Calculated MW	37
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only