

NKG2D Polyclonal Antibody Biotin Conjugated

Catalog No: #C90362Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NKG2D Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC-P IHC-F IF
Species Reactivity	Hu Ms Rat
Immunogen Type	121-219/219
Conjugates	Bio
Target Name	NKG2D
Other Names	NKG2D
Accession No.	Uniprot:O54709Gene ID:27007
Uniprot	O54709
GeneID	27007
Calculated MW	25
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only