

## ADRB2 Polyclonal Antibody Biotin Conjugated

Catalog No: #C90364Bio

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	ADRB2 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC P IHC-F IF
Species Reactivity	Hu Ms Rat
Immunogen Type	201-300/418
Conjugates	Bio
Target Name	ADRB2
Other Names	ADRB2
Accession No.	Uniprot:P07550Gene ID:154
Uniprot	P07550
GeneID	154
Calculated MW	46
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only