

phospho-NFKB p65 (Ser536) Polyclonal Antibody Cy3 Conjugated



Catalog No: #C90376Cy3

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	phospho-NFKB p65 (Ser536) Polyclonal Antibody Cy3 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	FS(p-S)IA
Conjugates	Cy3
Target Name	phospho-NFKB p65 (Ser536)
Other Names	phospho-NFKB p65 (Ser536)
Accession No.	Uniprot:Q04206Gene ID:5970
Uniprot	Q04206
GeneID	5970
Calculated MW	61
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only