

# CD71/Transferrin receptor Polyclonal Antibody PE-Cy5 Conjugated

Catalog No: #C90377PE-Cy5

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CD71/Transferrin receptor Polyclonal Antibody PE-Cy5 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt
Species Reactivity	Hu Ms Rat
Immunogen Type	621-720/760
Conjugates	PE-Cy5
Target Name	CD71/Transferrin receptor
Other Names	CD71/Transferrin receptor
Accession No.	Uniprot:P02786Gene ID:7037
Uniprot	P02786
GeneID	7037
Calculated MW	74/85
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

Flow-Cyt=1:50-200

Note: This product is for in vitro research use only