

CTLA4 Polyclonal Antibody PE Conjugated

Catalog No: #C90383PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CTLA4 Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	75-170/223
Conjugates	PE
Target Name	CTLA4
Other Names	CTLA4
Accession No.	Uniprot:P16410Gene ID:1493
Uniprot	P16410
GeneID	1493
Calculated MW	21
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=2ug/Test ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only