

EGFR Polyclonal Antibody AP Conjugated

Catalog No: #C90384AP

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	EGFR Polyclonal Antibody AP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC P IHC-F ICC
Species Reactivity	Hu Ms Rat
Immunogen Type	51-150/1210
Conjugates	AP
Target Name	EGFR
Other Names	EGFR
Accession No.	Uniprot:P00533Gene ID:1956
Uniprot	P00533
GeneID	1956
Calculated MW	130
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200

Note: This product is for in vitro research use only