

Occludin Polyclonal Antibody FITC Conjugated

Catalog No: #C90387FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Occludin Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	151-250/522
Conjugates	FITC
Target Name	Occludin
Other Names	Occludin
Accession No.	Uniprot:Q16625Gene ID:100506658
Uniprot	Q16625
GeneID	100506658
Calculated MW	59
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=2ug/Test ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only