

NCR1 Polyclonal Antibody Biotin Conjugated

Catalog No: #C90391Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NCR1 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB
Species Reactivity	Hu Ms Rat
Immunogen Type	198-206/304
Conjugates	Bio
Target Name	NCR1
Other Names	NCR1
Accession No.	Uniprot:O76036Gene ID:9437
Uniprot	O76036
GeneID	9437
Calculated MW	31
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200

Note: This product is for in vitro research use only