

Von Willebrand Factor Polyclonal Antibody FITC Conjugated



Catalog No: #C90398FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Von Willebrand Factor Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	1651-1800/2813
Conjugates	FITC
Target Name	Von Willebrand Factor
Other Names	Von Willebrand Factor
Accession No.	Uniprot:P04275Gene ID:7450
Uniprot	P04275
GeneID	7450
Calculated MW	309
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only