

ACE2 Polyclonal Antibody APC Conjugated

Catalog No: #C90399APC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ACE2 Polyclonal Antibody APC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	601-700/805
Conjugates	APC
Target Name	ACE2
Other Names	ACE2
Accession No.	Uniprot:Q9BYF1Gene ID:59272
Uniprot	Q9BYF1
GeneID	59272
Calculated MW	87
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only