

MX1 / MX2 Polyclonal Antibody Gold Conjugated

Catalog No: #C90405Gold

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MX1 / MX2 Polyclonal Antibody Gold Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IEM IGS GICA
Species Reactivity	Hu Ms Rat
Immunogen Type	301-400/662
Conjugates	Gold
Target Name	MX1 / MX2
Other Names	MX1 / MX2
Accession No.	Uniprot:P20591 Gene ID:4599
Uniprot	P20591
GeneID	4599
Calculated MW	73/79
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IEM=1:20-200 IGS=1:20-200 GICA=1:20-200

Note: This product is for in vitro research use only