

HGFA Inhibitor 2 Polyclonal Antibody PE Conjugated

Catalog No: #C90406PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	HGFA Inhibitor 2 Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu Ms Rat
Conjugates	PE
Target Name	HGFA Inhibitor 2
Other Names	HGFA Inhibitor 2
Accession No.	Uniprot:O43291Gene ID:10653
Uniprot	O43291
GeneID	10653
Calculated MW	28
Concentration	1 mg/ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only