

SFTPC Polyclonal Antibody APC Conjugated

Catalog No: #C90407APC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	SFTPC Polyclonal Antibody APC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rat
Immunogen Type	24-58/197
Conjugates	APC
Target Name	SFTPC
Other Names	SFTPC
Accession No.	Uniprot:P11686Gene ID:6440
Uniprot	P11686
GeneID	6440
Calculated MW	44307
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only