

# Dopamine Receptor D1 Polyclonal Antibody PE-Cy5.5 Conjugated



Catalog No: #C90414PE-Cy5.5

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Dopamine Receptor D1 Polyclonal Antibody PE-Cy5.5 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	101-200/446
Conjugates	PE-Cy5.5
Target Name	Dopamine Receptor D1
Other Names	Dopamine Receptor D1
Accession No.	Uniprot:P21728Gene ID:1812
Uniprot	P21728
GeneID	1812
Calculated MW	50
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only