

CDKN1A/P21 Polyclonal Antibody PE Conjugated

Catalog No: #C90424PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CDKN1A/P21 Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	21-100/164
Conjugates	PE
Target Name	CDKN1A/P21
Other Names	CDKN1A/P21
Accession No.	Uniprot:P38936Gene ID:1026
Uniprot	P38936
GeneID	1026
Calculated MW	18
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only