

Integrin alpha 2/CD49b Polyclonal Antibody Biotin Conjugated

Catalog No: #C90426Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Integrin alpha 2/CD49b Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB ICC
Species Reactivity	Hu Ms Rat
Immunogen Type	801-900/1181
Conjugates	Bio
Target Name	Integrin alpha 2/CD49b
Other Names	Integrin alpha 2/CD49b
Accession No.	Uniprot:P17301Gene ID:3673
Uniprot	P17301
GeneID	3673
Calculated MW	127
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 ICC=1:50-200

Note: This product is for in vitro research use only