

# Tetanus toxin light chain Polyclonal Antibody HRP Conjugated

Catalog No: #C90433HRP

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Tetanus toxin light chain Polyclonal Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC-P IHC-F ICC
Species Reactivity	Clostridium tetani E88
Immunogen Type	365-470/1315
Conjugates	HRP
Target Name	Tetanus toxin light chain
Other Names	Tetanus toxin light chain
Accession No.	Uniprot:P04958Gene ID:24255210
Uniprot	P04958
GeneID	24255210
Calculated MW	50
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200

Note: This product is for in vitro research use only