

FLT3L Polyclonal Antibody PE Conjugated

Catalog No: #C90437PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	FLT3L Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	61-160/235
Conjugates	PE
Target Name	FLT3L
Other Names	FLT3L
Accession No.	Uniprot:P49771Gene ID:2323
Uniprot	P49771
GeneID	2323
Calculated MW	23
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only