

FoxP3 Polyclonal Antibody HRP Conjugated

Catalog No: #C90447HRP

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	FoxP3 Polyclonal Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC-P IHC-F ICC
Species Reactivity	Hu Ms Rat
Immunogen Type	331-431/431
Conjugates	HRP
Target Name	FoxP3
Other Names	FoxP3
Accession No.	Uniprot:Q9BZS1Gene ID:50943
Uniprot	Q9BZS1
GeneID	50943
Calculated MW	47
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200

Note: This product is for in vitro research use only