

FoxP3 Polyclonal Antibody PE Conjugated

Catalog No: #C90447PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	FoxP3 Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	331-431/431
Conjugates	PE
Target Name	FoxP3
Other Names	FoxP3
Accession No.	Uniprot:Q9BZS1Gene ID:50943
Uniprot	Q9BZS1
GeneID	50943
Calculated MW	47
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=0.2ug/Test ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only