

CD5 Polyclonal Antibody PE Conjugated

Catalog No: #C90451PE



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	CD5 Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	221-320/495
Conjugates	PE
Target Name	CD5
Other Names	CD5
Accession No.	Uniprot:P06127Gene ID:921
Uniprot	P06127
GeneID	921
Calculated MW	55
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1ug/Test ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only